



Are the MicroSolv AQ brand syringe filters sterile – FAQ

No, but on a custom order we can supply some of our syringe filters in individual blister packs that can be sterilized.

Depending on the membrane, the sterilization could be autoclaving or radiation. Click [HERE](#) for ordering information and pictures of AQ™ syringe filters

Attachment

MicroSolv Filters Equivalency Study pdf [Download File](#)

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 09-26-2024