



Are RSA autosampler vials pyrogen free – FAQ

YES. The RSA™ or Reduced Surface **Activity** autosampler vials are pyrogen free when manufactured, as they are produced sterile and are never touched by “human hands”. They are then packaged, including shrink wrapped in a sterile, class 10,000 clean room environment.

We do not classify these vials as sterile for sale since the packaging is not sterilized before use. However, this can give you a very good idea of the cleanliness if you wished to sterilize them.

The entire manufacturing process is ISO Certified.

[More information about RSA vials, inserts and caps](#)



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 09-26-2024