

## Sample loop information – FAQ

**Sample Loops are available in a variety of internal volumes.**

Some of these will have several turns of the loop but others will not have a loop at all. In the latter case, the reason is because the tubing length is too small to make a loop. Whether or not there is a loop will not affect the functionality of the product.



Click [HERE](#) for sample loops ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)