



Autosampler vials sold by MicroSolv will fit properly in your instruments - Tech Information

All Autosampler vials sold by MicroSolv will fit commercially available modern HPLC autosamplers.

During the vial manufacturing processes, opto-electronic devices in the vial converting machines check within fractions of a second to ensure the processed vials meet the physical specifications (fit and form, dimension, etc.). In the case of a mismatch, the vial is blown out into the trash.

This 100 percent automatic control in-process controls as well as a final inspection according to DIN / ISO standards, which further ensures functionality and a perfect fit in the instrument. Worry free that the vials will cause missed injections or jam the autosamplers.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 09-19-2024