

---

## **Autosampler Vials sold by MicroSolv will Fit Your Instruments.**

All Autosampler Vials sold by MicroSolv will fit commercially available modern HPLC Autosamplers.

During the vial manufacturing processes, opto-electronic devices in the Vial converting machines check within fractions of a second to ensure the processed Vials meet the physical specifications (fit and form, dimension, etc.). In the case of a mismatch, the vial is blown out into the trash.

This 100 percent automatic control manual in-process controls as well as a final inspection according to DIN/ISO standards, which further ensures functionality and a perfect fit in the instrument.