



## Autosampler vials sold by MicroSolv will fit properly in your instruments – Tech Information

**All Autosampler vials sold by MicroSolv will fit commercially available modern HPLC autosamplers.**

During the vial manufacturing processes, opto-electronic devices in the vial converting machines check within fractions of a second to ensure the processed vials meet the physical specifications (fit and form, dimension, etc.). In the case of a mismatch, the vial is blown out into the trash.

This 100 percent automatic control in-process controls as well as a final inspection according to DIN / ISO standards, which further ensures functionality and a perfect fit in the instrument. Worry free that the vials will cause missed injections or jam the autosamplers.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 03-05-2024