



What is the purpose of a partial slit septa – FAQ

The partial slit septa offers both a penetration aid for thick and dull needles plus prevents evaporation before injection resulting in concentration changes. Only the silicone rubber is cut, the PTFE is not pre cut acting as an evaporation barrier.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 05-05-2024