



Wikipedia reference for use of pipets with organic solvents – Primer

From Wikipedia

Piston-driven air displacement pipettes

These are the most accurate and precise pipettes. They are more commonly used in biology, though they are commonly used by chemists as well.

The plastic pipette tips are designed for aqueous solutions, and are not recommended for use with organic solvents which may dissolve the plastic.

[Link to Wikipedia page](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com