



Are sparging stones good for ozonating water – FAQ

Yes, you can use any of our stainless steel sparging stones with ozone to “**ozonize**” water. It is recommended to use 10um pore size stones.

[*Click here for ordering information and pictures of our sparging stones*](#)

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 05-17-2024