



## Composition of linearity solutions in the HSQ Kit - Tech Information

The HSQ Kit's™ **linearity** solutions are manufactured under GMP conditions with calibrated balances and following a strict protocol and comprise of **caffeine** dissolved in 90:10 **acetonitrile**:0.1% **acetic acid**.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-22-2024