



Composition of linearity solutions in the HSQ Kit – Tech Information

The HSQ Kit's™ **linearity** solutions are manufactured under GMP conditions with calibrated balances and following a strict protocol and comprise of **caffeine** dissolved in 90:10 **acetonitrile**:0.1% **acetic acid**.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 10-11-2024