



Composition of linearity solutions in the HSQ Kit – Tech Information

The HSQ Kit's™ **linearity** solutions are manufactured under GMP conditions with calibrated balances and following a strict protocol and comprise of **caffeine** dissolved in 90:10 **acetonitrile**:0.1% **acetic acid**.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com