



How many injections should I expect on the HSQ Kit HPLC column - FAQ

Thousands... if you do not drop it, overheat it, put too much pressure through it over 4,000 psi, subject it to *pressure shock* or inject anything on the column that is not part of the HSQ Kit™.

Click [HERE](#) for HSQ Kit ordering information and pictures.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-22-2024