



I currently use a glass filter why should I use the HSQ Kit solutions instead for qualification - FAQ

Glass filters for UV detectors are not considered primary standards whereas the solutions in the MicroSolv HSQ Kit™ are primary standards.

Click [HERE](#) for PQ and HSQ ordering information and pictures.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-22-2024