



How many tests will I get out of one HSQ Kit for HPLC qualification - FAQ

Depending on how careful you are with transferring the solutions, each kit will produce **up to 30 tests**. These tests can be on one instrument or on 30 different instruments.

Click [HERE](#) for PQ Kit™ and HSQ Kit™ ordering information and pictures.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-22-2024