



How many tests will I get out of one HSQ Kit for HPLC qualification – FAQ

Depending on how careful you are with transferring the solutions, each kit will produce **up to 30 tests**. These tests can be on one instrument or on 30 different instruments.

Click [HERE](#) for PQ Kit™ and HSQ Kit™ ordering information and pictures.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com