

**Yes, this is why it can be run overnight from an autosampler.**

Simply set your conditions, program the system for performance qualification (PQ) validation and let it run. Use the following conditions:

Wavelength: 273nm  
Injection Volume: 10 $\mu$ l (except for volume **linearity**)  
Flow Rate: 1 mL/minute  
Temperature: Ambient

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)