



## If I am using TFA in HPLC will I have a problem with my Delrin type fittings - FAQ

**Delrin is not compatible with trifluoroacetic acid (TFA).**

If it comes in contact with your fittings, your fitting could become compromised and your connection might slip and become less than perfect or acceptable.

Click [HERE](#) to view our HPLC fittings

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)