



## Do all PEEK HPLC fittings fit all receiving ports – FAQ

No, you must consider and match up to the receiving port if it is Metric or English, the thread diameter, the “pilot length” of the receiving port, the thread count or pitch and the pressure expected.

*If you are not sure if your fitting should be used with a receiving port, contact the fitting and receiving port manufacturer for advice.*

Click [HERE](#) for HPLC fittings ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)