



Which internal diameters are the controlled flow CE capillaries available - FAQ

The Zero Flow™ and Controlled Flow™ capillaries are available in both 50m and 75m ID internal diameter. The length of the capillaries are meter.

Click [HERE](#) for Zero and Controlled Flow capillaries ordering information and pictures.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-03-2024