



## Which internal diameters are the controlled flow CE capillaries available – FAQ

The Zero Flow™ and Controlled Flow™ capillaries are available in both 50m and 75m ID internal diameter. The length of the capillaries are meter.

Click [HERE](#) for Zero and Controlled Flow capillaries ordering information and pictures.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 06-05-2024