



## Will all MicroSolv headspace crimp caps work for all brands of 20mm crimp top vials - FAQ

**Yes, most definitely.**

However, the crimp top (serum finish) for MicroSolv vials are tested for fit and form to insure precise crimping each and every time. Always be sure to have your crimping tool calibrated approximately every 50 crimps.

For more information about calibrating your crimping tool, please contact Technical Service for assistance.

[Click here for more information about Headspace Vials and Caps](#)

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)