



Which crimp cap is most commonly used for residual solvent testing by the pharmaceutical industry - FAQ

The crimp cap with silicone rubber / PTFE septa is most used by the pharmaceutical industry to be in compliance with USP <467>.

Our catalog number [95025-09-1S](#) is one of many caps used. Click [HERE](#) for more information about Headspace vials and caps

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com