



What is the advantage of using a Pharmafix type of septa for headspace analysis – FAQ

The advantage is that you get a better seal with the “Pharmafix” seal. However, you should be aware that it may not be as inert as a silicone/PTFE septa designed for headspace analysis.

Click [HERE](#) for ordering information and pictures of headspace caps

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com