



What is the purpose of foil as part of a septa in headspace analysis - FAQ

The purpose of the foil is generally to be a barrier film. To keep any solvents or gases from affecting the septa in the crimp cap. It also has a heat insulator property.

Click [HERE](#) for more information about Headspace Caps

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 06-16-2024