



## What is the purpose of foil as part of a septa in headspace analysis – FAQ

The purpose of the foil is generally to be a barrier film. To keep any solvents or gases from affecting the septa in the crimp cap. It also has a heat insulator property.

Click [HERE](#) for more information about Headspace Caps

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)