

What are Cogent HPS HPLC Columns – FAQ

Cogent HPS™ is a stationary phase (*both bonded and non bonded*) that is designed for applications that require high purity silica.

Metal content in Cogent HPS™ is extremely low and this produces very good stability and “deactivates” the silica for use with basic analytes. A very smooth and uniform silica surface results for mechanical strength and stability.

Cogent HPS™ is available for HPLC with unmodified silica, C18, C8, Amino, Cyano or Phenyl stationary phases.

[Click HERE for More Information about Cogent HPS HPLC Columns](#)

