

## What is the pH stability of Cogent HPS HPLC columns - FAQ

**These columns are extremely stable in from pH 2.0 to pH 8.6.**

Depending on the concentration of the acid or base will determine if these columns would be stable outside of this range.

*[Click here to view more information about Cogent HPS HPLC columns](#)*



Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-18-2024