

What are the metals content in Cogent HPS silica - FAQ

Metal	ppm Levels specification	Found
Na	<10	<1
Al	<10	<1
Fe	<10	<1
Ti	<10	<1
K	<10	None detected



Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-22-2024