



## What are the knitted reactor coils KRC manufactured with - FAQ

The knitted reactor coils KRC are made of virgin, medical grade PTFE (a very pure form of Teflon™).

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-22-2024