



## Are there any solvents that I cannot use with the knitted reactor coils KRC and knitted delay tubes KDT - FAQ

Yes, you cannot use any perfluorinated solvents or any solvent not compatible with PTFE with the knitted reactor coils KRC and knitted delay tubes KDT.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)