



Are there any solvents that I cannot use with the knitted reactor coils KRC and knitted delay tubes KDT - FAQ

Yes, you cannot use any perfluorinated solvents or any solvent not compatible with PTFE with the knitted reactor coils KRC and knitted delay tubes KDT.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 02-05-2024