



## What materials are the Cogent Column Filters manufactured with - FAQ

**The device (body) is made from high purity, inert PEEK with an internal stainless steel filter element and these are the only materials in the pre-column, Cogent Column Filters.**

*[Click here to view more information about pre-column filters for HPLC](#)*



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)