



What materials are the Cogent Column Filters manufactured with – FAQ

The device (body) is made from high purity, inert PEEK with an internal stainless steel filter element and these are the only materials in the pre-column, Cogent Column Filters.

[Click here to view more information about pre-column filters for HPLC](#)



Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 05-08-2024