



What is the maximum pressure the Cogent Column Filters should be used with – FAQ

The pre-column, Cogent Column Filters can be safely used with a maximum pressure of 5,000 psi.

[Click here to view more information about pre-column filters for HPLC](#)



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com