



## What is the maximum pressure the Cogent Column Filters should be used with – FAQ

**The pre-column, Cogent Column Filters can be safely used with a maximum pressure of 5,000 psi.**

*[Click here to view more information about pre-column filters for HPLC](#)*



Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 03-05-2024