



Can the Cogent Column Filters replace my HPLC guard columns - FAQ

No, a guard column is both a particulate filter (non dissolved chemicals) and a chemical filter (removes dissolved chemicals). The *Cogent Column Filter* is a particulate filter and adds extra protection for your expensive analytical and guard columns.

Click [HERE](#) for column protection ordering information.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com