

**No, a guard column is both a particulate filter (non dissolved chemicals) and a chemical filter (removes dissolved chemicals).** The *Cogent Column Filter* is a particulate filter and adds extra protection for your expensive analytical and guard columns.

Click [HERE](#) for column protection ordering information.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)