



Can the Cogent Column Filters replace my HPLC guard columns - FAQ

No, a guard column is both a particulate filter (non dissolved chemicals) and a chemical filter (removes dissolved chemicals). The *Cogent Column Filter* is a particulate filter and adds extra protection for your expensive analytical and guard columns

Click [HERE](#) for Column Protection Ordering Information



MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com