



What is the thickness of the SS filter element in the Cogent Column Filters - FAQ

The stainless steel filter element or “screen” is approximately 0.039” thick in the pre-column, Cogent Column Filters.

[Click here to view more information about pre-column filters for HPLC](#)



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com