



## Are the MicroSolv AQ brand syringe filters sterile – FAQ

**No, but on a custom order we can supply some of our syringe filters in individual blister packs that can be sterilized.**

Depending on the membrane, the sterilization could be autoclaving or radiation. Click [HERE](#) for ordering information and pictures of AQ™ syringe filters

### **Attachment**

**MicroSolv Filters Equivalency Study pdf** [Download File](#)

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-16-2024