



## Using one syringe filter for many samples will I get sample carry over - FAQ

**You will most likely get some kind of carry over.** It is not recommended to use one filter for more than one sample.

*[Click HERE for AQ syringe filter ordering information and pictures](#)*

### **Attachment**

**MicroSolv filters equivalency study pdf** [Download File](#)

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)