

Can I autoclave the MicroSolv AQ brand 25mm PTFE syringe filters - FAQ

Yes, under conditions of 121°C for a maximum of 15minutes, the PTFE AQ™ brand of syringe filters can be autoclaved. Maximum operating temperature (during continuous continuous autoclaving) we suggest that you do not exceed 80°C to maintain optimum performance.

Click [HERE](#) for ordering information and pictures of AQ™ brand of PTFE syringe filters

Attachment:

MicroSolv Filters Equivalency Study pdf [Download File](#)