



Are the AQ brand syringe filters bi-directional in their flow – FAQ

Yes, all the AQ™ Syringe Filters have bi-directional flow capabilities except the PTFE filter as well as other syringe filters that have glass fiber pre-filters.

Although in any syringe filter format, the flow in two directions may not be equivalent because of the filter device flow design. Click [HERE](#) for AQ brand syringe filter ordering information and pictures.

Attachment:

MicroSolv AQ filter equivalency study [Download File](#)

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 05-05-2024