



## Does the replaceable heating module come with the MicroSolv Window Maker when you order it or does it have to be ordered separately - FAQ

When you order the MicroSolv Window Maker™, it is shipped with a heating module that creates a detection window that is approximately 2mm. The replacement module is available in the event you need to replace it or other size window modules are available in 1mm, 4mm or 7mm.

Click [HERE](#) for Window Maker ordering information and pictures.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-21-2024