



Are there other parts or pieces that I need to use the MicroSolv Window Maker - FAQ

No, the MicroSolv Window Maker™ comes with a heating module that makes a 2mm window. No other tools are needed. Other heating modules of different window sizes can be purchased separately.

Click [HERE](#) for MicroSolv Window Maker ordering information and pictures.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-22-2024