



Are there other parts or pieces that I need to use the MicroSolv Window Maker - FAQ

No, the MicroSolv Window Maker™ comes with a heating module that makes a 2mm window. No other tools are needed. Other heating modules of different window sizes can be purchased separately.

Click [HERE](#) for MicroSolv Window Maker ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com