

Can I buy tubes made with RSA glass – FAQ

RSA™ glass is not available in any tube format and at this time.

RSA™ or Reduced Surface **Activity** glass is only available in 12x32mm autosampler vials. These vials will not release any contaminants such as chloroform or any other solvent into your sample.

For larger and smaller size vials, you may want to consider RSA-Pro X surface treatment. This is a non stick, hydrophobic surface treatment. Click [HERE](#) for information about RSA-Pro X vials and custom surface treatments.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com