



## Is The Vial Centrifuge NRTL rated – FAQ

No, not at this time. The Vial Centrifuge™ is not rated under the Nationally Recognized Testing Laboratories (NRTL) however, it is CE marked which indicates that the product complies with European health, safety, and environmental production standards.

Click [HERE](#) for ordering information.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-27-2024