



What is best way to know if there is a crack in a CE capillary - FAQ

The best way to know is to look at the entire capillary under either a stereo microscope or using a hand held magnifier offered on our website.

Click [HERE](#) for ordering information and pictures of the magnifier for this purpose.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 07-22-2024